

SMART PORTABILITY
CAN PROVIDE MORE VALUE






Get benefits with the best quality x-ray images

- Excellent image quality using direct deposited CsI
- Ultimate sharpness image by TRUVIEW™ART
- Instant upgrade to digital mobile x-ray system
- High resistance to impact and vibration
- Low price fixed grid (120 lines/inch)


1. High sharpness images by direct deposited CsI

Directly deposited CsI can provide clearer images at lowest dispersion compared to a conventional CsI and GOS scintillator. RadmediX promises its users best image quality with its design of low noise TFT, electronics and optimized CsI thickness.



2. Better DQE performance in higher spatial frequencies


RadmediX with its directly deposited columnar structured CsI + TFT provides outstanding high quality image with high DQE. It also demonstrates comparably excellent DQE performance in high spatial frequency range.



*DQE measurement condition : RQA-5 (Typical 2.6 uGy)


3. Lossless AED

Conventional AED function must consist of three steps: x-ray sensing, panel rest, and charge accumulation. Integration time is delayed as extra time is required for panel reset which occurs after the panel's x-ray sensing is achieved. Even when separate sensor modules within the detector system are used, x-ray loss is inevitable due to panel reset. With thick objects the loss rate can increase more. More accurate and stable AED function is possible with the Acuity DR Series. Lossless AED innovatively improved the reliability of sensitivity through operating scheme optimization.




4. Faster image display by a high speed operating scheme

Acuity DR is second to none in image preview and display time. Image preview in both wired and wireless conditions achieve times less than 2 seconds and image display is achieved in 4.6 and 6 seconds, respectively.




5. The best quality and x-ray image fine tuning with AccuVueMED & TRUVIEW™ART

AccuVueMED (X-ray UI software) and TRUVIEW™ART (image processing engine) use reverse filtering technology to decrease light dispersion, improve image sharpness in indirect type image processing where blurring occurs due to light dispersion. TRUVIEW™ART which comes standard in AccuVueMED enables more accurate diagnosis with high quality and high definition image.



Slim Dual Battery Charger

Two Acuity DR 1417 batteries can be simultaneously charged with the Slim Dual Battery Charger. LED display allows for easy checking of the battery level and with its light and compact design, user portability and usability is increased. The charger comes standard with the Acuity DR 1417 detector.



- Easy to handle with a battery knob design
- Battery check button and LED indicator


Durability

Mechanically optimized and engineered with metal unibody structure, the Acuity DR 1417 has safety at 3.3 ft height drop test. Also, protection suit is separately provided to protect the detector in various situations such as hospitals and mobile environments. Ergonomically designed protection suit provides optimum usability in x-ray environments. X-ray environment efficiency considered package of a protection suit and a compatible tablet mounter further adds to user convenience.



WPCS – Wireless Power Charging System

The Acuity DR 1417 is embedded with a wireless charging system that enables seamless 24 hours wireless operating environment. This was made possible with a detector embedded RX and external battery charging system TX. WPCS technology of the Acuity DR 1417 makes possible operation with wireless charging without the need of a battery change. WPCS is a world first new innovative product which installed a wireless battery charging system to extreme portability.



Inside AP & PCP

Inside AP maximizes the portability function of the Acuity DR 1417 which directly communicates with a Tablet PC and Smart Phones to make possible the viewing of images on mobile devices. With PCP (Portable Console PC) and wireless communication, the Acuity DR 1417's high definition images can be given anywhere and anytime. This enables outdoor and mobile applications such as in ICU, veterinary and industrial situations. Now, capture x-ray images on your own mobile devices.




PCP



Inside AP






CE 0123 UL Health Canada Class 2

Wifi & Wired Transfer

- Slim cassette size (43 x 43 cm, 15mm)
- Highly reliable and stable operation
- High definition by direct deposition CsI
- WiFi image transfer available (option)
- Fast image acquisition time less than 2 sec.
- Lossless AED support.



Excellent Post-processing Image Quality

Optimized algorithms and parameters for each body part
Adaptive noise reduction to minimize image signal loss
Image detail enhancement by multi-frequency image processing

User Experience Design





Smart workflow minimizing the need for page switches and mouse click
Editable tool bar
Easy to use stitching (up to 5 images)

Image Parameter Tuning Wizard

Users can select from 9 image styles processed using different parameters on a 3x3 matrix display.

Tablet, Smart Phone Supports

Supports viewing of high-definition digital images on display devices with WiFi communication such as Smart Phones and Tablets.





TRUVIEWTMART

ADVANCED IMAGE RECONSTRUCTION TECHNOLOGY





Image sharpness of an object in a conventional image is reduced due to light dispersion, TRUVIEWTMART, RadmediX's unique reverse filtering technology, reconstructs and improves image sharpness to increase the possibility of detecting abnormalities.

* Patent Pending





Specifications

| Model | Acuity DR 1417 | Acuity DR 1717 | | | | | | | | |
|-------------------------|--|--|-----------------------------|----------|---------------------------|--|-------|-----------------------------|----------|---------------------------|
| Detector Image |  |  | | | | | | | | |
| Purpose | Gen. Rad./ Wireless Portable | Gen. Rad. | | | | | | | | |
| Detector Type | CsI by direct deposition | CsI by direct deposition | | | | | | | | |
| Dimensions | 386 (H) x 460 (V) x 14.5 (D) | 460 (H) x 460 (V) x 15.0 (D) | | | | | | | | |
| Weight | 2.98 kg | 4.5 kg | | | | | | | | |
| Active Area | 358 x 430 mm | 430 x 430 mm | | | | | | | | |
| Pixel Pitch | 140 um | 140 um | | | | | | | | |
| Resolution | 2,560 x 3,072 | 3,072 x 3,072 | | | | | | | | |
| A/D Conversion | 14 bits | 14 bits | | | | | | | | |
| Image Aquisition | <table border="1"> <tr> <td>Wired</td> <td>Preview 2 sec/ Full 4.5 sec</td> </tr> <tr> <td>Wireless</td> <td>Preview 2 sec/ Full 6 sec</td> </tr> </table> | Wired | Preview 2 sec/ Full 4.5 sec | Wireless | Preview 2 sec/ Full 6 sec | <table border="1"> <tr> <td>Wired</td> <td>Preview 2 sec/ Full 4.5 sec</td> </tr> <tr> <td>Wireless</td> <td>Preview 2 sec/ Full 7 sec</td> </tr> </table> | Wired | Preview 2 sec/ Full 4.5 sec | Wireless | Preview 2 sec/ Full 7 sec |
| Wired | Preview 2 sec/ Full 4.5 sec | | | | | | | | | |
| Wireless | Preview 2 sec/ Full 6 sec | | | | | | | | | |
| Wired | Preview 2 sec/ Full 4.5 sec | | | | | | | | | |
| Wireless | Preview 2 sec/ Full 7 sec | | | | | | | | | |
| Voltage | DC 12V, 5 A | DC 12V, 5 A | | | | | | | | |
| Power Consumption | 12 W maximum | 12 W maximum | | | | | | | | |
| Communication Interface | 1 Giga bps/ IEEE 802.11n (2.4/5 GHz) | 1 Giga bps/ IEEE 802.11n (2.4/5 GHz) | | | | | | | | |
| X-ray Interface | Lossless AED | Lossless AED/ Wired Interface Cable | | | | | | | | |
| Accessory Parts | System Synchronization Unit Battery Charger (EVS-BCS) Battery Pack (EVS-MBP) | System Synchronization Unit | | | | | | | | |
| Certification |  |  | | | | | | | | |