




Acuity

FDR+



WORKFLOW INNOVATION

Radmedix 2510 Lance Drive, Kettering, OH 45409 USA
Tel: 844.723.6334 | Fax: 513.278.0150 | info@radmedix.com
© 2023 Radmedix Inc. All rights reserved.



FDR+ ADVANCED SYSTEM

Our Acuity FDR+ x-ray system with a standard configuration of a floor mounted tube stand, vertical bucky stand, and bucky table is shown below. This system features an advanced, 10" user-friendly touchscreen that provides a seamless operation of your radiographic system. It displays graphical tube angulation, SID values, and tube to bucky auto sync. With a range of advanced safety workflow and positioning features, the Acuity FDR+ is a high performance x-ray system offering unparalleled price-performance ratio.

10" Intuitive Touchscreen

Preset source to image distance with Auto Sync feature.

Real time display of SID and tube angulation.

Suitable for any clinical application.



Premium Wall Stand

Motorized and counterbalance adjustment controls with electromagnetic locks.

Convenient handgrips to optimize patient positioning.

Rotational bucky to easily change the panel orientation.



AEC Optional

FDR+ KEY HIGHLIGHTS



Auto Sync Tracking with Wall & Table Bucky's




Rotational Column



6 Way Elevating Float Top Table




Rotational Bucky



AEC Optional

Emergency Stop




Wireless Foot Switch

EFFORTLESS POSITIONING

Motorized auto-tracking enables a wide array of patient positioning.



INTUITIVE 10" TOUCHSCREEN



Elevated Table Specifications

- Lowering Height: 21 in.
- Patient Load: 550 lbs
- Tabletop dimensions: 32 in. x 86 1/2 in.
- Integrated Emergency Stop

FDR STANDARD SYSTEM

The Acuity FDR is also offered in our standard configuration. Offering the same versatility of the Acuity FDR+, but with our standard non-touchscreen display and advanced features. The Acuity FDR Standard offers the same exceptional digital imaging quality but with standard features at a budget friendly price point.

4" Console Display

Preset source to image distance with Auto Sync feature.

Real time display of SID and tube angulation.

Suitable for any clinical application.



AEC Optional

Premium Wall Stand

Manual counterbalance adjustment controls with electromagnetic locks.

Convenient handgrips to optimize patient positioning

Rotational bucky is only available on the FDR+.



FDR STANDARD FEATURES



Rotational Column



4 Way Float Non Elevating Table



Built In Foot Pedal



AEC Optional